



Filter Reliability

All Semrock filters demonstrate exceptional reliability. The simple all-glass structure and hard dielectric glass coatings (as hard as the glass on which they are coated!) mean they are virtually impervious to humidity-induced shifting and can be cleaned and handled like any standard glass optics. All Semrock coatings are hard – **never** soft.

All Semrock filters are capable of withstanding high optical intensities. Semrock filters either have laser-damage threshold specifications, or, depending on the application, have been tested with intense broadband sources for extended periods of time with no noticeable degradation. **Semrock filters do not burn out** under normal conditions, even with prolonged use.


You will never find any adhesives in the optical path of a Semrock filter. Epoxies and cements can absorb water vapor and swell or shift, and they can photo-darken or be optically damaged. **Semrock filters are built for longevity.**

Semrock filters have been tested to meet or exceed requirements for environmental and physical durability set forth in specifications such as MIL-STD-810F, MIL-C-48497A, MIL-C-675C, and ISO 9022-2. The table below shows some of the key standards against which our filters are regularly tested.

Environmental Durability Testing	Standard/Procedure	Test Description
Humidity	MIL-STD-810F (507.4)	Aggravated Humidity (> 10 x 24 hr cycles)
High Temperature	MIL-STD-810F (501.4)	Induced Hot (> 7 x 24 hr cycles)
Low Temperature	MIL-STD-810F (502.4)	Cold (C2) (24 hr cycles)
Physical Durability Testing	Standard/Procedure	Test Description
Adhesion	MIL-C-48497A (4.5.3.1)	"Tape Test"
Humidity	MIL-C-48497A (4.5.3.2)	Damp Heat
Moderate Abrasion	MIL-C-48497A (4.5.3.3)	"Cheesecloth Test" (> 50 cycles)
Solubility/Cleanability	MIL-C-48497A (4.5.4.2)	Sequential immersion in acetone and alcohol
Water Solubility	MIL-C-48497A (4.5.5.3)	Immersion in distilled water (> 24 hrs)



Semrock's dedicated reliability laboratory with state-of-the-art environmental test chambers is in regular use.



Examples of bandpass filter spectra measured before and after ten 24-hour cycles of Aggravated Humidity testing according to MIL-STD-810F, demonstrating negligible change.