




### Technical Information: RazorEdge® Filter Layouts

Filters play a crucial role in enabling high-fidelity Raman spectroscopy measurements. In many simple laboratory layouts the only filters needed are a [laser-transmitting filter](#) to clean-up the laser spectral output and a [laser-blocking filter](#) to ensure that no Rayleigh-scattered laser light reaches the detection system (see [Filter Types for Raman Spectroscopy Applications](#)). However, in many high-performance Raman systems – such as those with microscopic imaging capabilities or highly sensitive remote probes – expensive and/or complex focusing and collection optics are used to couple the system to the sample region. For these systems it is desirable for the excitation laser beam and the Raman-shifted signal light to share a common light path. The illustrations below show how this layout can be accomplished with an ultrasteep dichroic beamsplitter used in conjunction with a laser-blocking filter. Semrock's [RazorEdge Dichroic™](#) beamsplitter filters are ideal for these applications.



*In the standard Raman spectroscopy layout, the laser excites the sample directly, and a laser-blocking filter (such as the ultrasteep RazorEdge filter) is positioned between the sample and the spectrometer with light incident at or near normal incidence – its job is to block the scattered laser light and pass the Raman-shifted signal light as close as possible to the laser wavelength and with as much transmission as possible. RazorEdge filters are ideal for this configuration.*

*In focusing or imaging systems that utilize high-NA collection optics, for example, it is convenient for the incident laser beam and collected signal light to share a common path. To meet this requirement, a two-filter solution is ideal: a 45° beamsplitter reflects the laser light and directs it through the optics to the sample, while efficiently transmitting the returning Raman-shifted signal light; a laser-blocking filter at normal incidence is used in conjunction with the beamsplitter to completely block the undesired laser light from the spectrometer.*

In order for the two-filter configuration shown above right to work, the 45° beamsplitter must be as steep as the laser-blocking filter. Traditionally thin-film filters could not achieve very steep edges at 45° because of the "polarization splitting" problem – the edge position tends to be different for different polarizations of light. However, through continued innovation in thin-film filter technology, Semrock has been able to achieve ultrasteep 45° beamsplitters with the same steepness of our renowned RazorEdge laser-blocking filters: the transition from the laser line to the passband of the filter is guaranteed to be less than 1% of the laser wavelength (for U-grade filters)!